

HV for hodoscope

output pulse for beta source (electrons pathing through) set to 200 ± 20 mV

channel	HV	channel	HV
1	1795 V	17	1475 V
2	1750 V	18	1775 V
3	1600 V	19	1575 V
4	1750 V	20	1675 V
5	1795 V	21	1550 V
6	1675 V	22	1400 V
7	1550 V	23	1725 V
8	1795 V	24	1725 V
9	1675 V	25	1675 V
10	1795 V	26	1700 V
11	1675 V	27	1666 V
12	1650 V	28	1795 V
13	1566 V	29	1575 V
14	1550 V	30	1475 V
15	1625 V	31	1650 V
16	1625 V	32	1600 V

red - lower amplitude (cca 2x)